

4L/4R CURETTE

128-002



PRODUCT DESCRIPTION

Double ended instrument - both ends a mirror image of one another to give left and right curvature

Features

- A double ended instrument - the two ends are a mirror image to give a left and right (complimentary) curvature.
- Curettes are modified hand scalers which have smaller toes (working area) with rounded tips.
- A balanced grip handle which has a large diameter for comfortable grip and good control.

Uses

- To scale subgingivally (remove calculus below the gumline and in periodontal pockets) and to perform root planning (remove redundant, necrotic cementum).
- They may also be used supragingivally to remove light calculus deposits
- Periodontal pockets (deeper than 2mm) should be cleaned with hand instruments (A curette) rather than power scalers.

How to use

- The instrument is held in modified pen grip
- The toe is placed against the tooth surface with the flat face at 90 degrees so that the sharp edge contacts.
- Repeat overlapping strokes are used.
- The movement is from the base of the periodontal pocket towards the tip of the crown.

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	128-002

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights