

**EXPLORER PROBE - PATTERN NO.6**

128-059

**PRODUCT DESCRIPTION**

**A sharp pointed probe with a robust yet flexible tip and high-quality handle for a comfortable grip**

**Key Features**

- A sharp pointed probe
- No. 6 pattern = simple, almost 90 degree, bend in the shank with a straight tip
- The tip is fine, slightly flexible yet robust.
- A deluxe quality handle gives a comfortable grip, good control and tactile response

**Used to feel for:**

- Defects in the tooth surface
- Pulp exposure sites in fractured teeth tooth decay (early cavities)
- Early feline resorptive lesions.
- Detect calculus subgingivally

**How to use:**

The instrument is held in a modified pen grip. Run the sharp point over the tooth surface.

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	128-059
<b>Generic Size</b>	X-Small

Colour/Style

Stainless Steel

Dental Instrument Type

Other

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights