

LEAK TESTERS FOR FLEXIBLE ENDOSCOPES

END-LEAK



PRODUCT DESCRIPTION

Leak testers are required for all immersible flexible endoscopes and leak testing should be performed prior to use and prior to immersion

Select the right leak tester for your endoscope, or call us for advice on cleaning and disinfecting your endoscopes

ADDITIONAL INFORMATION

SKU	END-LEAK
Animal	Mixed

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

