

## PASSIVE SCAVENGING KIT

103-114



### PRODUCT DESCRIPTION

## Easy to install passive scavenging kit requiring no power

A simple method of removing the waste anaesthetic gases, directing the gases outside the building with an easy to install passive scavenging wall outlet

- The outlet point should be installed approximately 400mm from ground level
- The tubing from the circuit to the outlet should not exceed 4 metres in length
- Optional exterior cawling available to protect the outside of the installation
- Includes internal stainless steel fascia plate for interior
- Installation available
- 22mm Passive gas scavenging elbow

### ADDITIONAL INFORMATION

**Animal**

Mixed

**SKU**

103-114

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

