

## SET OF 5 X LOW DEADSPACE CONNECTORS

VSL-VS0055



### PRODUCT DESCRIPTION

---

## Reduce dead space and rebreathing with your sidestream capnography connectors

Although often required, attachments used between a patient and a breathing system can lengthen the apparatus dead space and contribute to the rebreathing of expired gases, especially when the fresh gas flow is not adequate. Heat and moisture exchangers can also increase this.

Our low dead space adapters allow a Sidestream capnography line to be directly attached, reducing dead space.

Set of 5 x RSP 15mm ISO ET Connectors with Sideport PK (2.6,3,3.6,5,6mm)

**Compatible with Mindray and Vetronic Monitors**

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	VSL-VS0055
<b>Type</b>	Other
<b>Monitor Fit</b>	Vetronic Impact, Vetronic Lightning
<b>Accessory Parameter</b>	CO <sup>2</sup>

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights