

## GIMA DIATERMO MB ELECTROCAUTERY UNIT & ACCESSORIES

GIMA-MB120D-BUNDLE\*



### PRODUCT DESCRIPTION

## A high-frequency electrosurgical machine for performing middle monopolar or bipolar electrosurgery.

Allowing pure cut (CUT), cut-coagulation (BLEND), superficial coagulation (FORCED COAGULATION), forced coagulation in the absence of necrosis (SOFT COAG) and, with forceps, Bipolar coagulation.

The GIMA MB 120 guarantees absolute reliability with monitoring of operative function and a digital reading of delivered power. The neutral electrode connection and patient/plate contact are constantly monitored throughout operation.

The possibility of function and power control via the handle reduces distractions and improves focus on the procedure.

#### Key Features:

- Minimum invasive surgery
- Compatible with bipolar forceps
- Memorises the last used settings.
- Power activation by a foot pedal or handle
- Digital adjustment and an indication of power output
- Adjustment of the acoustic indicator intensity.

#### Technical Specifications:

- Minimum preselect able power: 0
- Level step: 1
- Maximum output power CUT (W): 120 250 Ohm
- Maximum output power BLEND (W): 90 200 Ohm
- Maximum output power COAG FORCED (W): 80 150 Ohm
- Maximum output power COAG SOFT (W): 60 100 Ohm
- Maximum output power BIPOLAR (W): 40 100 Ohm
- Working frequency: 600 kHz
- Weight kg: 7
- Size LxHxD mm: 260 x 110 x 265
- Operating voltage (switchable): 115 -230

- Mains frequency (Hz): 50-60
- Fuses 230Vac (5x20) TIMED: 2 x 3,15A
- Fuses 115Vac (5x20) TIMED: 2 x 6,3A
- Electrical input power (VA): 300
- Electrical input current (230 Vac) (A): 1,3
- Electrical input current (115 Vac) (A): 2,6
- Five steps adjustable sound level (from 55- to 75dBA): YES
- Self-check: YES
- Power accuracy output warning: YES
- Split or not split patient plate allowed: YES
- Last working condition storing: YES
- Electrical Class (EN60601 - 1): I CF
- MDD 93/42/EEC Class: II B
- Patient circuit: F

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	GIMA-MB120D-bundle*
<b>Cut Modes</b>	pure cut, coagulated cut, forced coagulation, soft coagulation, bipolar coagulation

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights