

GIMA DIATERMO MB 120D MONO-BIPOLAR ELECTROSURGERY UNIT - CLEARANCE

GIMA-MB120D-CL



PRODUCT DESCRIPTION

A HF electrosurgical unit for performing middle monopolar or bipolar electrosurgery.

Includes FREE Autoclavable MB Handle (GIMA-30554)

The modes of operation they provide are pure cut, coagulated cut, forced coagulation, soft coagulation and, with forceps, bipolar coagulation.

The use of either split of normal neutral plate is allowed.

Leading technology enables highly professional surgery intervention thanks to ergonomic, plus safety features derived from high power theatre apparatus.

An open circuit alarm warns the user that connection into the neutral plate has not been made or that the cable is disconnected or the neutral plate is not in contact with the patient. The alarm is signalled by an intermittent beep, together with a flashing red light on the front panel.

OVT alarm: the unit includes an internal circuit protection warning when the current delivery is exceeding.

Technical Specifications:

- Minimum preselect able power: 0
- Level step: 1

https://burtonsveterinary.com/ Page 1 Sales: 01622 834300 | Service: 01622 834350



- Maximum output power CUT (W): 120 250 Ohm
- Maximum output power BLEND (W): 90 200 Ohm
- Maximum output power COAG FORCED (W): 80 150 Ohm
- Maximum output power COAG SOFT (W): 60 100 Ohm
- Maximum output power BIPOLAR (W): 40 100 Ohm
- · Working frequency: 600 kHz
- Weight kg: 7
- Size LxHxD mm: 260 x 110 x 265
- Operating voltage (switchable): 115 -230
- Mains frequency (Hz): 50-60
- Fuses 230Vac (5x20) TIMED: 2 x 3,15A
 Fuses 115Vac (5x20) TIMED: 2 x 6,3A
- Electrical input power (VA): 300
- Electrical input current (230 Vac) (A): 1,3
- Electrical input current (115 Vac) (A): 2,6
- Five steps adjustable sound level (from 55- to 75dBA): YES
- · Self-check: YES
- Power accuracy output warning: YES
- Split or not split patient plate allowed: YES
- · Last working condition storing: YES
- Electrical Class (EN60601 1): I CF
- MDD 93/42/EEC Class: II B
- · Patient circuit: F

S/N: 0269009

ADDITIONAL INFORMATION

Animal	Mixed
SKU	GIMA-MB120D-cl
Cut Modes	pure cut, coagulated cut, forced coagulation, soft coagulation, bipolar coagulation

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights