

## BURTONS CARBONTECH TABLE MAT V2

165-434-TABLEMAT



### PRODUCT DESCRIPTION

---

## Our CarbonTech heat mats feature next-generation technology for maximum efficiency and safety.

Our latest design, the V2 includes 4 extra safety features, a safe temperature regulation sensor, 3 selectable power levels, shutdown after 12 hours of use, and a disconnectable power block. This table size mat measures to 112cm L x 52cm W.

The carbon technology uses high-density carbon fabric to evenly distribute the heat across the entire surface of that mat. Its soft material will move easily into any shape, providing maximum patient comfort and coverage.

Once switched on, it will reach its target temperature in less than 5 minutes and will maintain that temperature whilst transferring the warmth direct to the patient. For this to happen, the patient must be in direct contact with the sensor. This sensor also works as a safe temperature regulator, not allowing the transferred temperature to exceed 42C.

CarbonTech heat mats are designed for use without any other materials between the patient and the blanket, having any other insulation, such as blankets or bedding can interfere with the temperature sensors, and therefore interfere with the correct temperature being transferred to the patient.

The mat does not store heat, so will significantly reduce the risks of hot spots and burns compared to other warming mats on the market when used in the correct way. The protective coating also provides maximum resistance to animal scratches and bites, and a 3-year warranty is provided for blanket material damage. It is also completely waterproof and easy to clean.

#### Key Features:

**Carbon Technology** - CarbonTech uses high-density carbon fabric specifically designed for even heat distribution.

**Flexible material** - The soft nature of the heat mat allows you to bend it around the patient, providing maximum patient comfort and coverage.

**Easy to use** - Once switched on, the mat will reach its target temperature in less than 5 minutes. CarbonTech is designed to then maintain its target temperature.

**Safe heating** - CarbonTech works by delivering a low-contact temperature whilst transferring only the energy necessary to warm the patient. It does not store heat, reducing hot spots.

**Safe Temperature Sensor** - Integrated sensor means the mat will only work when the patient makes direct contact with it, it also contains a limiter, not allowing the temperature transfer to exceed 42C.

**Disconnect-able power bank** - This updated feature means less disturbance for your patients. If you want to leave your patient on the mat once its not in use to save disturbing them, you can simply disconnect the entire power cable and bank, allowing you to get rid of any hanging wires and use the mat without its heating element.

**3 Power Levels** - Set which power level you need before turning the heat mat on.

**Robust** - The protective coating on CarbonTech is designed for maximum resistance to animal scratching and biting, we provide a 3-year warranty on patient damage to the mat.

**Hygienic** - Can be cleaned easily with water and soap. After drying, it can be disinfected with alcohol.

**Safety shutdown** - The heat mat will automatically turn itself off after 12 hours, if it is for any reason, accidentally left on.

## ADDITIONAL INFORMATION

|               |                  |
|---------------|------------------|
| <b>SKU</b>    | 165-434-TABLEMAT |
| <b>Brand</b>  | Mano Medical     |
| <b>Animal</b> | Small Animal     |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights