

## MEDICAPTURE 300 VIDEO ARCHIVING UNIT (HD)

END-MCUSB300



### PRODUCT DESCRIPTION

---

**The new MediCap USB300 offers several new enhancements, building upon its reputation as the world's highest quality and easiest to use medical video recorder**

New features include an ULTRA recording mode with twice the quality of Blu-Ray, 3D recording capability and 1TB storage. Making it the perfect choice for high performance recording.

#### Key Features:

- **HD/3D Medical Video Recorder** - Records HD and 3D video from your endoscope or surgical camera. Saves to an internal hard drive, external USB drive, or network drive.
- **Intuitive Use** - It is so easy, no training is needed. Just plug it in and go. You will be recording videos within minutes of opening the box.
- **Medical Grade** - With its steel case and fluid resistant front panel, the USB300 is built to withstand the rigors of the OR, and it is so compact it will fit on any cart

#### Specifications:

- **Device Type:** HD/3D2 Medical Video Recorder Records high-definition video and still images from medical video sources to internal hard drive, external USB drive, or network drive
- **Case:** Rugged metal case with high-impact plastic faceplate
- **LCD Monitor:** 2 1/2", full colour (for video confirmation only, not for diagnosis)
- **Buttons:** Sealed membrane, fluid resistant
- **Video Inputs:** DVI (RGB via adapter), HD SDI, S-Video, Composite

- **Video Formats:** H.264 / MPEG4 (1080, 720, NTSC, PAL)
- **Video Quality:** Compact, Normal, Ultra (up to 40 Mbit/s)
- **Image Formats:** JPG, TIFF, BMP, DICOM
- **Networking:** RJ45 10/100/1000 Ethernet
- **Internal Storage:** Internal hard drive, 1TB
- **External Storage:** USB flash drive, USB hard drive, network drive
- **USB Support:** USB2.0, high speed, 1 front connector, 2 rear connectors
- **Networking:** Ethernet T10/100 connection via Network Upgrade Kit (optional)
- **Footswitch/ Trigger:** Hands-free capture Standard 3.5mm mini jack
- **Mechanical:** Size: 240 x 210 x 63mm / 9.5 x 8.27 x 2.5" Weight: 1.7kg / 3lb 12oz (appx.)
- **Environmental:** Storage: -40° to +85° C Operating: -20° to +40° C
- **Power:** AC 100-240V 50/60Hz 3-prong AC jack
- **Media:** USB flash drives or external USB hard drive
- **Image/ video resolution:** NTSC: 720 x 480 pixels PAL: 720 x 576 pixels
- **Image Formats:** JPEG, TIFF, BMP, DICOM
- **Video format:** H.264 (DVD-quality)
- **Audio Format:** AAC, 2 channels
- **Certifications:** CE 93/42/EEC, IEC 60601-1:2005, ANSI/AAMI ES60601-1:2005, CAN/CSA C22.2 No.60601-1:2008, E361228, FCC Part 15, FDA Class 1, RoHS, PSE, GOST-R
- **Kit Contents:** MediCap USB200, USB flash drive, video cables, power cable, user guide

## ADDITIONAL INFORMATION

<b>SKU</b>	END-MCUSB300
<b>Animal</b>	Mixed
<b>Generic Size</b>	Small
<b>Colour/Style</b>	White

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights