

## ENDOSCOPY STORAGE CABINET

END-ESCGB



### PRODUCT DESCRIPTION

**Store your endoscopes securely with this lockable endoscope storage cabinet**

**Key Features:**

- Endoscope storage cabinet holds 2 scopes
- Bracket included
- Dimensions 2150mm H x 500mm W
- For bespoke cabinets call our sales team

### ADDITIONAL INFORMATION

|               |           |
|---------------|-----------|
| <b>SKU</b>    | END-ESCGB |
| <b>Brand</b>  | Burtons   |
| <b>Animal</b> | Mixed     |

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights