


**OVAL LOOP GRABBER SNARE 1.8MM X 240CM 2.5CM LOOP**

END-F026

**PRODUCT DESCRIPTION**

**Our Oval Loop Grabber Snare has inner barbed wire for effective retrieval in flexible procedures**

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	END-F026
<b>Channel Size</b>	1.8mm
<b>Autoclavable</b>	Yes
<b>Working Length</b>	2400mm

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

