

## END-ICAPTSC TOUCHSCREEN COMPUTER WITH SOFTWARE

END-ICAPTSC



### PRODUCT DESCRIPTION

---

**I-CAP touchscreen computer preloaded with i-CAP V4 and required operating systems for still and video image capture**

#### Key Features

- Touchscreen PC with keyboard and mouse
- Export images to included report templates or create your own
- Capture, catalogue and report
- Designed for endoscopy and MIS
- 1 data entry to capture images
- Ideal for PPE in clinic, hospitals and universities

#### Extent Of Delivery

- Main Unit
- Keyboard
- Mouse
- Mains power cable
- Manual
- Video cables
- Installation

### ADDITIONAL INFORMATION

**SKU**

END-ICAPTSC

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights