

## NON-AUTOCCLAVABLE WATER BOTTLE FOR OLYMPUS ENDOSCOPES

END-WB-140



### PRODUCT DESCRIPTION

Generic replacement water for the 140, 160, 180 Series and newer Olympus Endoscopes.

This replacement water bottle has a 250ml capacity and offers the following features:

- Easy to connect/disconnect from scope
- Metal connector to bottle top and metal endoscope connector for durability
- Attached cap to protect endoscope connector
- Metal clip and strong plastic strap to attach the water bottle to the stock hanger on Olympus light sources
- Non-autoclavable - use soaking methods to sterilise

### ADDITIONAL INFORMATION

SKU

END-WB-140

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

