

## REPLACEMENT BLOOD PRESSURE CUFFS FOR PETMAP

129-561-1



### PRODUCT DESCRIPTION

## Replacement blood pressure cuffs for the PetMap range of monitors

Criter cuffs are veterinary blood pressure (BP) cuffs specifically designed for use with petMAP range of monitors. The markings found on the cuffs allow for easy selection of the correct size for each animal. Made with high-quality materials, Criter cuffs are very durable.

- Available in range of sizes from 2.0cm - 13.0cm
- Extension available

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	129-561-1
<b>Colour/Style</b>	White
<b>Brand</b>	PetMAP
<b>Catagory Sub Heading</b>	Range of Sizes Available
<b>Type</b>	Cuff
<b>Monitor Fit</b>	PetMAP Graphic

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights