

## IMEC DISPOSABLE TEMPERATURE EXTENSION CABLE

IMEC-0011-30-37391



### PRODUCT DESCRIPTION

## Temperature extension cable for disposable probes

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	IMEC-0011-30-37391
<b>Brand</b>	Mindray
<b>Type</b>	Lead
<b>Monitor Fit</b>	Mindray ePM12, Mindray iMEC 8, Mindray iPM12, Mindray uMEC12
<b>Accessory Parameter</b>	Temperature

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

