

**MONITOR SAMPLE LINE**

V8044

**PRODUCT DESCRIPTION****8ft (2.4m) CO<sup>2</sup> Sample Line**

- Priced per line
- Fit: V8200, V8401, V8402, V90041, V90043, V9400, Advisor monitor and most other side stream CO<sub>2</sub> monitors with a luer connector
- Only available to order in quantities of 10

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	V8044
<b>Brand</b>	SurgiVet
<b>Type</b>	Other
<b>Monitor Fit</b>	Surgi Vet Advisor

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights