

SAMPLE LINE - 8FT(2.4M) CO²

V8044

**PRODUCT DESCRIPTION****8ft (2.4m) CO² Sample Line**

Fits: V8200, V8401, V8402, V90041, V90043, V9400, SurgiVet Advisor monitor and most other side stream CO² monitors with a luer connector

Change every 2-4 weeks. Compatible with M880B, Edan X12, Mindray uMEC, Mindray ePM12 & SurgiVet Advisor.

ADDITIONAL INFORMATION

| | |
|------------------------------|--|
| Animal | N/A |
| SKU | V8044 |
| Accessory Parameter | CO ² |
| Monitor Compatibility | BLT M880B, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor |
| Change Frequency | 2-4 weeks |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights