

NASAL SAMPLE LINE - 10 PK

V1186



PRODUCT DESCRIPTION

Nasal sample line for capnographs

- For monitors: V8200, V8401, V8402, V990041, V90043, V9400 and most other side stream CO2 monitors with a luer connector.

ADDITIONAL INFORMATION

Animal	N/A
SKU	V1186
Brand	SurgiVet
Type	Other
Monitor Fit	Burtons Capnograph, Burtons M880, Mindray iMEC 8, Mindray iPM12, Surgi Vet Advisor

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

