

SURGIVET RE-USEABLE PURPLE CUFF

31543B-CONFIG



PRODUCT DESCRIPTION

Comfortable fit design

These re-useable Surgivet cuffs are design to fit patients comfortably thanks to its 360° bladder that eliminates pinching and also eliminates the need to align an arrow/mark over an artery.

The durable cuffs can withstand 10,000 inflations to 200 mmHg and are machine washable. A reinforced closure system ensures that the cuffs will stay on and will handle the wear and tear of everyday practice use.

- Suitable for SurgiVet V6004 and Advisor Monitors

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	31543B-config
Colour/Style	Purple
Brand	SurgiVet
Catagory Sub Heading	Range of Sizes Available
Type	Cuff
Monitor Fit	Surgi Vet Advisor

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights