

## SURGIVET RE-USEABLE PURPLE CUFF

31543B-CONFIG



### PRODUCT DESCRIPTION

## Comfortable fit design

These re-useable Surgivet cuffs are design to fit patients comfortably thanks to its 360° bladder that eliminates pinching and also eliminates the need to align an arrow/mark over an artery.

The durable cuffs can withstand 10,000 inflations to 200 mmHg and are machine washable. A reinforced closure system ensures that the cuffs will stay on and will handle the wear and tear of everyday practice use.

- Suitable for SurgiVet V6004 and Advisor Monitors

### ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	31543B-config
<b>Colour/Style</b>	Purple
<b>Brand</b>	SurgiVet
<b>Catagory Sub Heading</b>	Range of Sizes Available
<b>Type</b>	Cuff
<b>Monitor Fit</b>	Surgi Vet Advisor

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights