

C-FUSOR PRESSURE INFUSER

V6403



PRODUCT DESCRIPTION

C-Fusor 500 mL pressure infusor

The crystal clear material allows for maximum visibility of fluid levels. It features an accurate and easy to read pressure gauge, quick, secure fastening and a large squeeze bulb for rapid and reliable pressure control.

Features

- Secure fastening
- Durable and easy to clean
- Allows for even pressure around the fluid container
- Three-way stopcock prevents air leaks

ADDITIONAL INFORMATION

Animal	Mixed
SKU	V6403
Brand	SurgiVet
Type	Cuff, Other
Monitor Fit	Surgi Vet Advisor

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights