

## C-FUSOR PRESSURE INFUSER

V6403



### PRODUCT DESCRIPTION

## C-Fusor 500 mL pressure infusor

The crystal clear material allows for maximum visibility of fluid levels. It features an accurate and easy to read pressure gauge, quick, secure fastening and a large squeeze bulb for rapid and reliable pressure control.

#### Features

- Secure fastening
- Durable and easy to clean
- Allows for even pressure around the fluid container
- Three-way stopcock prevents air leaks

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	V6403
<b>Brand</b>	SurgiVet
<b>Type</b>	Cuff, Other
<b>Monitor Fit</b>	Surgi Vet Advisor

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights