

**ECG PRE-GELLED ELECTRODE PADS**

115-964

**PRODUCT DESCRIPTION****Ecg pre-gelled electrodes**

- Paediatric cloth backing
- Suitable for short term monitoring, less than 24 hours
- Pack of 3
- 19X38mm
- Solid gel, MRI/radiolucent

**ADDITIONAL INFORMATION**

<b>Animal</b>	N/A
<b>SKU</b>	115-964
<b>Brand</b>	SurgiVet
<b>Type</b>	Other
<b>Monitor Fit</b>	Other
<b>Accessory Parameter</b>	ECG

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights