

## HINGED BLADE CASE

AND-12370



## PRODUCT DESCRIPTION

- Keeps your blades safe and in good condition while on the move as a mobile groomer or in the salon
- The Andis Hinged Blade Case is completely portable
- Holds up to 12 blades
- Maximum protection with high-density foam

## BLADES NOT INCLUDED

## ADDITIONAL INFORMATION

SKU

AND-12370

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights