

## ECG TRUNK CABLE, 3/5 LEAD, CLIP, ADULT / PAEDIATRIC 12 PIN

IMEC-0010-30-42719



### PRODUCT DESCRIPTION

## ECG Trunk Lead

- 3/5 Lead
- Adult or pediatric use
- 12 Pin
- Defibrillator proof
- AHA / IEC

, 12 Pin For Mindray uMEC, Mindray ePM12 monitors.

### ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	IMEC-0010-30-42719
<b>Brand</b>	Mindray
<b>Accessory Parameter</b>	ECG
<b>Monitor Compatibility</b>	BLT E30VET, Edan X12, Mindray ePM12, Mindray uMEC, SurgiVet Advisor
<b>Change Frequency</b>	When Broken

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights