

STANDARD 2 TAP BATH/SHOWER MIXER IN CHROME

115-302



PRODUCT DESCRIPTION

Swan standard 2 tap with shower mixer for tub - chrome

Designed with a focus on functionality with high quality, hard wearing features, backed up by 5 year guarantee

- Hot & cold mixing and spout diverter
- Chrome finish
- Includes 1.5m hose, shower head and wall bracker for shower
- Minimum supply pressure 0.2 bar
- Material: brass

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	115-302
Operation	Hand

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

