

BURTONS 'CLASSIC' TWIN DEPTH TUB TABLE

125-205



PRODUCT DESCRIPTION

Our classic style twin depth tub table

The classic style split tub with 17.25-inch (43.81cm) deep end and 5.25-inch (13cm) shallow tub end in an over knee design. Giving you greater versatility for treating a variety of patients.

Overall length 60-inch (152cm) width 24-inch (61) height 36 inch (91cm). Tap fittings and drain sold separately.

Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality materials this product has been engineered from the ground up to last!

Features

- Supplied with a removable recessed grill
- Shower tap and drain fittings sold separately
- Overall Dimensions: 153cm x 61cm x 90cm LWH
- Tub Depths: 13cm/44cm

ADDITIONAL INFORMATION

| | |
|---------------|--------------|
| Animal | Small Animal |
| SKU | 125-205 |

Colour/Style

Stainless Steel

Brand

Burtons

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights