

BURTONS PRO-LIFT MONO TABLE

225-207-BUNDLE



PRODUCT DESCRIPTION

The latest addition to our range of tables, this ultra-quiet, electric examination table excels in functionality, stability and longevity

Get your patient to the correct position with ease with the Burtons Pro-Lift Mono. With adjustable height, controlled using the electric powered ram, reduce repetitive strain and ensure maximum comfort. Its simple controls enable you to adjust the table height from 690mm to 1090mm and the built-in manual tilt allows you to get optimum angle and positioning for all examinations.

With a standardised table top fitting this table can be configured for your needs allowing for most of the Burtons table top range to be fitted with additional heated and lead-lined options available. Simply select the top you need.

Electric Height Control

Enjoy ease of use around your practice with built-in battery, allowing you to work wherever you need, and don't be afraid of it losing power, the table can be used directly from mains as a backup!

The Frame

Made up of a single, heavy duty, electric ram the simple design of this table base allows for quick and easy cleaning. All supported

on a heavy duty stainless steel base with four anti-static medical grade, lockable castors, ensuring smooth and easy manoeuvrability around the practice.

Choose your top:

Standard Top

Standard flat top constructed of high-quality stainless steel. Welded and polished corners allow the top to be thoroughly wiped clean with no hidden dirt traps. Complete with tie down points to secure the patient. Top Size: 122 x 60 cm

Extended Top

Extended version of our flat top allowing additional patient / working space on the table. The perimeter of the top has tie down points to secure the patient. Top Size: 137 x 60 cm

V Top

For procedures requiring precise animal positioning, choose a 'V' top style table top. Featuring two die formed stainless steel table leaves, which can be positioned independently to help you create the required angle for specific procedures. Top Size: 147 x 40 cm

Dental Top

Perfectly designed for any dental procedure, the high-quality, stainless steel top is intelligently shaped to ensure any excess water runs into the built-in grill. Designed at an extended top size allowing for the dental top to be fitted to a range of Burtons table bases and also providing ample procedure space. Top Size: 137 x 60 cm

Heated Top - Max temp: 35°

The heated top option available on a range of table top assists in maintaining body temperature during use. The temperature is controlled via an easy-to-use mounted panel, simply select the maximum temperature you want the table top to reach and the sensor within the table will then regulate the temperature once reached - 137cm x 60cm

Lead Lining

Available on a range of tops, lead-lined tops are a high-quality stainless steel top with built-in lead lining, an ideal combination for X-ray procedures and operations/examinations where space is at a premium in the practice.

Don't forget to protect your top with a Burtons table mat..

We use a 2.5mm thick black matting made from a manmade material that is non-marking and has excellent chemical and abrasive resistance. Mats are provided with stainless steel securing clips for the two short ends, these help hold the mat in place on tables without a recess for the mat to lay in.

Manufactured to last

We don't just supply it! It is designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent. Using the latest CNC's and laser manufacturing processes and high-quality stainless steel materials this product has been engineered from the ground up to last!

Key features:

- Electric Lift
- Hand Controls
- Battery As standard
- Manual Tilt (10° increments)

- Hygienic and easy to clean
- Ultra-quiet operation
- Easy-to-use controls
- Heavy-duty, high-quality stainless steel base.
- Twin wheels 6" lockable, anti-static castors.
- Minimum height: 69cm
- Maximum height: 109cm
- Maximum load: 200kg
- Various flat and speciality top options

Dimensions:

- Electric Pro-Lift Table Base (844 x 400 x 1049) LxWxH

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	225-207-bundle
Colour/Style	Stainless Steel
Brand	Burtons
Catagory Sub Heading	Build Your Perfect Table
TiltingTop	Yes - Hydraulic
Mobile	Yes
Control Type	Handheld
Lift Type	Electric

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights