


BURTONS LO-LIFT TABLE

125-150-BUNDLE



Stable at every height



Locking system with 3 options;
Full lock, straight wheel only, free movement

PRODUCT DESCRIPTION

With its smartly designed, heavy duty lift mechanism, this table top will go just 47cm from the floor, allowing patients to be easily placed on the table and safely lifted to an optimum working height

The Lo-Lift table series built by Burtons has been designed around the Vet, with ultra-low minimum height, 3 pedal braking, electric battery option and high quality frame. This table is a staple of many veterinary clinics. With a standardised table top fitting this table can be configured for your needs allowing for most of the Burtons table top range to be fitted with additional heated and lead lined options available. Simply select the top you need.

Electric or Hydraulic

Choose from either the simple pedal pump hydraulic lift option or electric lift with built in battery and hand control system making even the heaviest dog easy to lift to the perfect working height.

The frame

The key part of this table is the superbly engineered, powder coated fully welded steel frame. Offering you stability throughout the

entire Lo-Lift range and huge load capacity of 180kg! Large, anti-static, heavy duty wheels go on to ensure smooth, safe mobility around the practice.

Three pedal brake for optimum manoeuvrability

When the pedal is totally depressed, the wheels cannot be moved at all, when the pedal is up, one wheel is fixed straight making the table easy to steer, and when the pedal is in a horizontal position, all the wheels are free.

Choose your top:

Standard Top -

Standard flat top constructed of high quality stainless steel. Welded and polished corners allow the top to be thoroughly wiped clean with no hidden dirt traps. Complete with tie down points to secure the patient. Top Size: 122 x 60 cm

Extended Top -

Extended version of our flat top allowing additional patient / working space on the table. The perimeter of the top has tie down points to secure the patient. Top Size: 137 x 60 cm

Heated Top - Max temp: 35°

The heated top option available on a range of table top assists in maintaining body temperature during use. The temperature is controlled via an easy to use mounted panel, simply select the maximum temperature you want the table top to reach and the sensor within the table will then regulate the temperature once reached

Lead Lining -

Available on a range of tops, lead lined tops are A high quality stainless steel top with built in lead lining, an Ideal combination for X-ray procedures and Operations/ Examinations where space is at a premium in the practice.

Don't forget to protect your top with a Burtons table mat

We use a 2.5mm thick black matting made from a manmade material that is non-marking and has excellent chemical and abrasive resistance. Mats are provided with stainless steel securing clips for the two short ends, these help hold the mat in place on tables without a recess for the mat to lay in.

Manufactured to last

We don't just supply it! It is designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

Key Features:

- Minimum height: 47cm
- Maximum height: 95cm
- Maximum load: 180kg
- Top: Various Flat top options

- Extra low height adjustment
- Lift: Hydraulic or Electric
- Battery: As standard with electric lift
- Controls: Hand controls (electric)
- Fully mobile with easy steer
- Unique wheel locking system
- Simple 4-wheel lock/brake
- Ideal all-rounder

Dimensions:

- Electric Lo-Lift (1189.5 x 699 x 798.1) LxWxH
- Hydraulic Lo-Lift (1189.5 x 699 x 798.6) LXWxH

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	125-150-bundle
Colour/Style	Cream
Brand	Burtons
TiltingTop	No
Mobile	Yes
Control Type	Footcontrol, Handheld
Lift Type	Electric, Hydraulic

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights