

## BURTONS OPERATING TABLE

225-099-BUNDLE



### PRODUCT DESCRIPTION

---

**With an extremely stable design, even when fully extended, this superbly built, stainless steel operating table with its easy clean boxed design, is the go-to operating table for many veterinary clinics.**

Adjusting the patient height to suit your needs is easy with either the hydraulic lift or electric lift options. Secure, quick release position locking levers enable you to alter the table top height and angle from horizontal to almost vertical and even rotate the top through 360° (Not available on electric lift option).

The table base is fitted with integrated 'wheel barrow' style wheels for superb stability, whilst allowing easy maneuvering whilst cleaning. With a standardised table top fitting, this table can be configured for your needs allowing for most of the Burtons table top range to be fitted with additional heated and lead-lined options available. Simply select the top you need.

#### Electric or Hydraulic Lift

Choose from either the simple pedal pump hydraulic lift option or electric lift with either hand or foot control system making even the heavier dog easy to lift to the perfect working height

## The frame

The key part of this table is the tubular ram in a fully welded, hygienic, stainless steel box design. Giving the table extra stability, 360° rotation off the top (hydraulic only) and easy cleaning.

## Wheelbarrow style maneuverability

Built with a two-wheel design meaning movement of the table only when you want it!

## 360° table top rotation

On Hydraulic Operating Table Base ONLY. Due to the internal workings of the electric operating table base, 360° rotation is not possible

## Choose your top:

### Standard Top

Standard flat top constructed of high-quality stainless steel. Welded and polished corners allow the top to be thoroughly wiped clean with no hidden dirt traps. Complete with tie-down points to secure the patient. Top Size: 122 x 60 cm

### Extended Top

An extended version of our flat top allowing additional patient/working space on the table. The perimeter of the top has tie-down points to secure the patient. Top Size: 137 X 60 cm

### V Top

For procedures requiring precise animal positioning, choose a 'V' style table. Featuring two die formed stainless steel table leaves, which can be positioned independently to help you create the required angle for specific procedures. Top Size: 147 x 40 cm

### Dental Top

Perfectly designed for any dental procedure, the high quality, stainless steel top is intelligently shaped to ensure any excess water runs into the built-in grill. Designed at a standard top size allowing for the dental top to be fitted to a range of Burtons table bases and also providing ample procedure space. Top Size: 137 x 60 cm

### Heated Top - Max temp: 35°

The heated top option available on a range of tabletop assists in maintaining body temperature during use. The temperature is controlled via an easy to use mounted panel, simply select the maximum temperature you want the table top to reach and the sensor within the table will then regulate the temperature once reached - 137cm x 60cm

### Lead Lining

Available on a range of tops, lead-lined tops are A high-quality stainless steel top with built-in lead lining, an Ideal combination for X-ray procedures and Operations/ Examinations where space is at a premium in the practice.

### Don't forget to protect your top with a Burtons table mat...

We use a 2.5mm thick black matting made from a manmade material that is non-marking and has excellent chemical and abrasive resistance. Mats are provided with stainless steel securing clips for the two short ends, these help hold the mat in place on tables without a recess for the mat to lay in.

## Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high-quality stainless steel materials this product has been engineered from the ground up to last!

## Key features

- Electric or hydraulic Lift
- Quick release levers, electric or foot controlled
- Hygienic and easy to clean
- Ultra-quiet operation
- High-quality stainless steel base.
- Minimum height: 80cm
- Maximum height: 111.5cm (including top)
- Maximum Tilt: 80°
- Maximum load: 136Kg
- Top - Various flat and specialty top options

## Dimensions:

- Hydraulic Operating Table Base (665 x 405 x 1070) LxWxH

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	225-099-bundle
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Catagory Sub Heading</b>	Build Your Perfect Table
<b>TiltingTop</b>	Yes - Hydraulic
<b>Mobile</b>	Yes
<b>Control Type</b>	Footcontrol, Handheld
<b>Lift Type</b>	Electric, Hydrualic

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights