

## BURTONS PRO-LIFT TWIN COLUMN TABLE

225-205-BUNDLE



### PRODUCT DESCRIPTION

---

## Never struggle to achieve the perfect patient position again with the newly designed, fully electric twin column Pro-Lift table

With a high capacity 400kg weight limit this ultra-quiet, electric examination table excels in functionality, stability, and longevity. Get even your heavier patient to the perfect position with ease using the fully electric height and tilt controlled auto-leveling dual rams. Reducing risk of repetitive strain and ensuring maximum comfort. Its simple controls enable you to adjust the table height from 690mm to 1090mm and tabletop tilt to 33°.

With a standardised table top fitting, this table can be configured for your needs allowing for most of the Burtons table top range to be fitted with additional heated and lead-lined options available. Simply select the top you need.

#### Fully electric base

Enjoy the ease of use around your practice with a built-in battery, allowing you to work wherever you need, and don't be afraid of it losing power, the table can be used directly from mains as a backup.

#### The frame

Made up of dual heavy duty, electric rams, the simple design of this table base allows for quick and easy cleaning. All supported on a heavy duty stainless steel base with four anti-static medical grade, lockable castors, ensuring smooth and easy manoeuvrability around the practice.

### **Choose your top:**

#### **Standard Top**

Standard flat top constructed of high quality stainless steel. Welded and polished corners allow the top to be thoroughly wiped clean with no hidden dirt traps. Complete with tie down points to secure the patient. Top Size: 122 x 60 cm

#### **Extended Top**

Extended version of our flat top allowing additional patient / working space on the table. The perimeter of the top has tie down points to secure the patient. Top Size: 137 X 60 cm

#### **V Top**

For procedures requiring precise animal positioning, choose a 'V' style table. Featuring two die formed stainless steel table leaves, which can be positioned independently to help you create the required angle for specific procedures. Top Size: 147 x 40 cm

#### **Dental Top**

Perfectly designed for any dental procedure, the high quality, stainless steel top is intelligently shaped to ensure any excess water runs into the built in grill. Designed at a standard top size allowing for the dental top to be fitted to a range of Burtons table bases and also providing ample procedure space. Top Size: 137 x 60 cm

#### **Heated Top - Max temp: 35°**

The heated top option available on a range of table top assists in maintaining body temperature during use. The temperature is controlled via an easy to use mounted panel, simply select the maximum temperature you want the table top to reach and the sensor within the table will then regulate the temperature once reached - 137cm x 60cm

#### **Lead Lining**

Available on a range of tops, lead lined tops are A high quality stainless steel top with built in lead lining, an Ideal combination for X-ray procedures and Operations/ Examinations where space is at a premium in the practice.

#### **Don't forget to protect your top with a Burtons table mat..**

We use a 2.5mm thick black matting made from a manmade material that is non-marking and has excellent chemical and abrasive resistance. Mats are provided with stainless steel securing clips for the two short ends, these help hold the mat in place on tables without a recess for the mat to lay in.

#### **Manufactured to last**

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

#### **Key Features:**

- Electric Lift and Tilt
- Hand Controls
- Battery as Standard
- Auto-levelling
- Hygienic and easy to clean

- Ultra-quiet operation
- Easy to use controls
- Heavy-duty, high-quality stainless steel base.
- Twin wheels - lockable, anti-static castors.
- Minimum height: 69cm
- Maximum height: 109cm
- Maximum Tilt: 33°
- Maximum load: 400Kg
- Top - Various flat and speciality top options

**Dimensions:**

- Pro-Lift Twin Table Base (899mm x 412mm x 1071mm) LxWxH

**ADDITIONAL INFORMATION**

<b>Animal</b>	Small Animal
<b>SKU</b>	225-205-bundle
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Catagory Sub Heading</b>	Build Your Perfect Table
<b>TiltingTop</b>	Yes - Electric
<b>Mobile</b>	Yes
<b>Control Type</b>	Handheld
<b>Lift Type</b>	Electric

**Good to know...**

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights