


DENTAL TOP FOR OPERATING TABLE

225-103-I



PRODUCT DESCRIPTION

The newly designed top has built in raised edges to ensure no water can run off the table when in use

Perfectly designed for any dental procedure, the high quality, stainless steel top is intelligently shaped to ensure any excess water runs into the built in grill.

Designed at a standard top size allowing for the dental top to be fitted to a range of Burtons table bases and also providing ample procedure space.

Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent. Using the latest CNC's and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

Top Specification

- Top Size: 137 x 60 cm

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	225-103-I

Colour/Style

Stainless Steel

Brand

Burtons

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights