

## 'V' TOP OPERATING TABLE

225-102\*VTOPCONFIG



### PRODUCT DESCRIPTION

---

**Featuring two die formed stainless steel table leaves, which can be positioned independently to help you create the required angle for specific procedures**

Features include Twist locking points at either end of the table for added stability and security and a built-in central runner tray under the center of the table allowing for simple control of any liquids to a drain/collection hole. Keeping the operating site clear. The table leaves can be easily removed for cleaning and a satin finish ensures this top is easy on the eyes even in bright theatre lighting!

Designed to be fitted to a range of Burtons table bases and also providing ample procedure space, this table top is ideally suited to our tilting range of operating tables, giving you effective control over the position of the patient.

- Independent leaf angle
- 6 tie down cleats as standard
- Built with a removable channel runner makes for easy cleaning and hygiene maintenance
- Top size 147 x 40 cm

#### **Why not make your top Heated ?**

The heated top option available on a range of tabletop assists in maintaining body temperature during use. The temperature is controlled via an easy to use mounted panel, simply select the maximum temperature you want the table top to reach and the

sensor within the table will then regulate the temperature once reached

### Choose your perfect base...

#### Burtons Operating Table base (Recommended)

Designed for maximum stability and easy cleaning the Burtons tilting operating base with different top options is perfect for operating procedures. With electric (hand or foot controlled) or hydraulic lift and 'wheelbarrow' style maneuvering.

#### Burtons Pro-lift Mono Table base

Configurable to fit your needs, the Pro-lift mono base with a configurable top is perfect for all types of examinations. It's high quality, ultra-quiet electric lift (up to 200kg) and manual 33° tilt makes it easily positionable and stable at any height!

#### Burtons Pro-lift Table Base

High quality fully electric lift and 33° electric tilt table base. The Pro-lift with its ultra-quiet, ultra smooth twin rams, and configurable tops is both stable and powerful with a maximum weight of 400kg! Perfect for all types of examinations and operations.

### Manufactured to last

We don't just supply it! It was designed and built by Burtons in the UK at our principal manufacturing facility in Maidstone, Kent Using the latest CNC's and laser manufacturing processes and high-quality stainless steel materials this product has been engineered from the ground up to last!

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	225-102*vtopconfig
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Catagory Sub Heading</b>	Build Your Perfect Table
<b>TiltingTop</b>	Yes - Hydraulic
<b>Mobile</b>	Optional
<b>Control Type</b>	Footcontrol, Handheld
<b>Lift Type</b>	Electric, Hydrualic

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights