

BURTONS CONSULTATION TABLE

115-023-BUNDLE



PRODUCT DESCRIPTION

***Due to international shipping disruption, black rubber matting is temporarily out of stock. We expect new stock mid-to-late next month.**

Orders are currently being fulfilled with white matting. For consulting tables, a free black mat will be sent when stock arrives if preferred. Replacement and cut-to-size orders will be held until black stock arrives unless you are happy with white.

Please [contact us](#) if you have any questions.*

A simple consult table that will last a lifetime, retain its looks, and impress your clients

Burtons consult tables are available in standard rigid or mobile versions. Supplied with a black table mat as standard, which contrasts smartly with the stainless steel frame and provides a warm, noise-free surface for the patient.

Finished in a brushed satin the table top is formed from a single sheet of stainless steel with welded and polished corners; this allows the top to be thoroughly wiped clean with no hidden dirt traps or rusting edges.

Key Features:

- 6 Tie down cleats available (optional extra)
- Dimensions: 116 x 55 x 84.15 (including static feet)
- Lead lined option
- Mobile or static base
- Optional Scales

- Scales - Weight range: 0 - 145kg (Optional extra)
- Scales - Increments: 50g (Optional extra)
- 150Kg max. static load

Configure your perfect consulting table

Keep you and your patients safe

6 x Securing cleats can be provided at a small additional cost if required.

Make it mobile

Mobile versions are fitted with 75mm diameter total lock swivel castors, which not only brake the table travel, but also hold it firmly in place, a particularly important feature when dealing with nervous animals.

Why not

Make your consultation table a set of scales for your practice by incorporating a weighing scale platform into the top. Ideal for consulting and prep rooms where space is at a premium.

The scale can be operated via the mains adaptor or from the internal rechargeable battery and has a useful weight hold feature included on the display

Save space and Combine your tables

Lead-lined versions are available as a further option, an Ideal combination for X-ray procedures and operations/examinations where space is at a premium in the practice.

- Scales - Weigh units: lbs or kg (Optional extra)

Manufactured to last

We dont just supply it! It is designed and built by Burtons in the UK at our principle manufacturing facility in Maidstone, Kent Using the latest CNCs and laser manufacturing processes and high quality stainless steel materials this product has been engineered from the ground up to last!

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	115-023-bundle
Colour/Style	Stainless Steel
TiltingTop	No
Mobile	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights