

## MOBILE CONSULTATION TABLE

115-025



## PRODUCT DESCRIPTION

A simple consult table that will last a lifetime, retain its looks, and impress your clients! Burtons' consulting tables are available in standard rigid or mobile versions. Supplied with a black tablemat as standard, which contrasts smartly with the stainless steel frame and provides a warm, noise-free surface for the patient. Finished in a brushed satin the table top is formed a single sheet of stainless steel with welded and polished corners; this means that the top can be thoroughly wiped clean with no hidden dirt traps or rusting edges.

Key Features:

- 6 Tie down cleats available (optional extra)
- Dimensions: 116x55x84.37 high (including castors)
- Lead lined option
- Mobile or static base
- Optional Scales
  - 150Kg max. static load

## ADDITIONAL INFORMATION

<b>Animal</b>	Small Animal
<b>SKU</b>	115-025
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>TiltingTop</b>	No
<b>Mobile</b>	No

## Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights