

BURTONS SCALE CONSULTING TABLE

115-047



PRODUCT DESCRIPTION

Ideal for consulting and prep rooms where space is at a premium

Make your consultation table a set of scales for your practice by incorporating a large weighing scale platform into the top. Ideal for consulting and prep rooms where space is at a premium.

The scale is ideal to weigh anything from cats to large dogs up to 150kg. The scale features Zero, Hold and Tare functions. Animals can be weighed in either kg/g or lb/oz.

Burtons' scales consulting tables are supplied with a black tablemat as standard, which contrasts smartly with the stainless steel frame and provides a warm, noise-free surface for the patient.

Key Features:

- Powered by batteries for use anywhere in the surgery
- Stainless Steel Platform
- Easy Clean, Hygienic Top
- Adjustable Feet
- Stainless steel frame

Specifications:

- Dimensions: 116x55x86cm
- Weight range: 0-145kg

- Increments: 50g
- Weigh Units : lbs or kg
- 1050mm x 505mm Scale Area
- 150Kg max. static load
-

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	115-047
Colour/Style	Stainless Steel
TiltingTop	No
Mobile	No

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.
Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!
Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.
Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights