

### 3 NAIL GRASPER - 1.8MM

END-10-3NG



## PRODUCT DESCRIPTION

3 Nail Graspers for foreign body retrieval - in various sizes

## ADDITIONAL INFORMATION

<b>Animal</b>	N/A
<b>SKU</b>	end-10-3NG
<b>Category Sub Heading</b>	Range of Working Lengths Available
<b>Channel Size</b>	1.8mm
<b>Autoclaveable</b>	Yes

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.  
**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!  
**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.  
**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

