

ELONGATED RAT TOOTH GRASPERS - 2.3MM

END-ELONGATED-RATTOOTH-2.3

**PRODUCT DESCRIPTION****Elongated rat tooth graspers for foreign body retrieval in flexible endoscopy procedures**

- Autoclavable
- 2 different lengths available
- For use in 2.3mm channels or larger

ADDITIONAL INFORMATION

Animal	Mixed
SKU	END-ELONGATED-RATTOOTH-2.3
Catagory Sub Heading	Range of Working Lengths Available
Channel Size	2.3mm
Autoclaveable	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights