

## 3 LOOP ATRAUMATIC GRASPERS - 2.3MM

END-LOOP-AT-2.3



### PRODUCT DESCRIPTION

## 3 Loop atraumatic graspers for gastroscopy procedures

- Autoclavable
- Ideal for foreign body retrieval
- For use in 2.3mm channels or larger
- Various working lengths available

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	END-LOOP-AT-2.3
<b>Category Sub Heading</b>	Range of Working Lengths Available
<b>Channel Size</b>	2.3mm
<b>Autoclavable</b>	Yes

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

