

## NSK HIGH SPEED HANDPIECE

148-550



### PRODUCT DESCRIPTION

**High-speed handpiece with a patented clean head system, designed to prevent oral fluid or foreign particles from entering into the handpiece head.**

The system is activated instantaneously in response to air pressure change as the handpiece stops spinning.

#### Key Features:

- Ceramic Bearings - long lasting
- Clean Head System for enhanced infection control
- Ultra Push Chuck - simple yet secure fastening
- Single Spray
- Speed : 380,000~480,000 min-1
- Head Size: ø10.3 x H13.8 mm
- Midwest 4 hole
- Single Spray - efficient cooling
- Ergo Grip - secure grip for your confidence
- Handpieces are autoclavable to 135°.

### ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	148-550
<b>Generic Size</b>	X-Small
<b>Examination Type</b>	Oral Examination

**Fitting Type**

4 Hole

**Burr Fit**

F/G

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights