

CHEEK DILATORS

128-DILATORS



PRODUCT DESCRIPTION

Cheek dilators for small furrries in 2 different sizes

Key Features

- Spring loaded – quick and safe application.
- Solid paddles – for added tissue protection.
- Stainless steel construction.

Uses:

- Small dilator - For small rabbits (under approx.1kg body weight), chinchillas and other small rodents.
- Large dilator - Used in combination with a mouth gag to retract the cheeks and so aid visualization of the mouth and improve access to the cheek teeth.
- Used in combination with a mouth gag
- To retract the cheeks and so aid visualization of the mouth and improve access to the cheek teeth
- To protect the inside of the cheeks during treatment, however other tissue guards need to be used in conjunction e.g wooden lollipop sticks

ADDITIONAL INFORMATION

Animal	Small Animal
SKU	128-dilators
Dental Instrument Type	Dilators

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights