

## WALL MOUNTED VETAIR WITH OPTIONAL COMPRESSOR

118-780-BUNDLE



### PRODUCT DESCRIPTION

---

**Want a VETair dental machine for your practice but struggling for floor space? The VETair wall-mounted system is the ideal solution when space is limited.**

The unit has a double arm to reach where you need it, but this also means that it folds back on itself minimizing the wall space it occupies. A removable tray is provided on the wall bracket. As with the mobile VETair this unit is supplied with curly leads so they don't drape across the floor and get damaged.

The water supply to the handpieces can also be used to supply your dental scaler, eliminating the need for additional water bottles and saving even more space.

The system requires an air supply. This is typically from a compressor (optional extra) which can be sited away from the point of use. As an alternative, the system can be adapted to use a medical air supply.

#### **Installation:**

It is important that the wall you intend to mount the system to is sufficiently strong to cope with the weight of the system. Wall fixings are not included as this will vary dependant on the wall structure. With each system we include a 6mm x 10mtr length of air supply tubing, longer lengths can be supplied on request. A separate air control box is included which needs to be sited within approx. 1m radius of the wall bracket, this could be hidden away, fixed in a cupboard or just wall mounted. The installation is very straightforward, a competent handyman will have no trouble installing the system. We will, of course, provide a commissioning visit to make sure everything is in order and to instruct users on the safe use and maintenance of the system.

**Why not add the VETair starter kit and one of our dental scalers to create a complete dental package.**

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	118-780-bundle
<b>Colour/Style</b>	Stainless Steel
<b>Brand</b>	Burtons
<b>Examination Type</b>	Oral Examination, Small Furries, Tooth Extraction

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights