

## HOSPIVAC 400 SUCTION UNIT

129-610



### PRODUCT DESCRIPTION

---

## The HOSPIVAC 400 has been designed for professional aspiration of bodily fluids, tissues or bones of the patient during or after surgery

Equipped with two 2000cc polycarbonate autoclavable jars with an overflow valve, the hospivac surgical aspirator is built with a state-of-the-art oilless pump guaranteeing high performances with excellent suction capacities and max vacuum reachable within a few seconds, with no need for maintenance. A clear dashboard along with a full range of accessories and ability to be continually used makes it the ideal device for longer surgical procedures.

The system is portable, mounted on four antistatic castors and features a transport handle so can be easily moved around the practice.

Supplied with suction tubing and disposable flexible suction end with optional foot pedal and rigid suction ends sold separately.

- Liquid Collection Jars with overflow valve system
- 8x14 Silicone Tubes, length 150 cm (autoclavable) + Conical Connector  $\varnothing$  10-11-12 mm
- Antibacterial & Hydrophobic Filter (single-patient)
- Power Cord
- 2x2000 ml Jars

### Technical Specifications

- Max Vacuum (adjustable) -0.90 bar -90 kPa -675 mmHg
- Max free air flow rate 90 l/min
- Noise Level 46,4 dBA
- Duty cycle Non-stop operation
- Weight kg 20
- Size cm 46 x 85 x 42

- Voltage - 230V / 50 Hz
- Power - 300 VA
- Duty Cycle - Continuous
- Motor Oilless and maintenance-free piston pump
- Complies with - Safety requirements for electrically powered suction equipment UNI EN ISO 10079-1 Class IIa Medical Device according to 93/42/EEC "Medical Device" Directive

## ADDITIONAL INFORMATION

<b>Animal</b>	Mixed
<b>SKU</b>	129-610
<b>Overall Product Dimensions (mm)</b>	470 x 740 (h) x 430
<b>Colour/Style</b>	White

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights