

## ASKIR SUCTION UNIT & ACCESSORIES

129-600-ACCESSORIES



### PRODUCT DESCRIPTION

## A compact high suction system ideal for the majority of veterinary applications including use with an endoscope

The unit is supplied with a 2-litre autoclavable collection jar which features a float valve to prevent overfilling. An anti-bacterial filter is fitted as standard and a 1.5-metre length of suction tubing is also supplied. The level of suction is controlled by a rotary valve and indicated by a pressure manometer. Suction ends sold separately. Disposable Suction End features a control port to limit suction. Place a finger over the port to increase suction.

#### Technical Information:

- Suction type: High vacuum/high flow
- Max Suction: -0.8 bar (-12psi)
- Max Free Air Flow rate: 40L/min
- Duty cycle: 120mins on / 60 mins off
- Jar capacity: 2 litres
- Dimensions: 37cmH x 21cmW x 35cmL
- Weight: 2.5kg
- Power: 230v/50Hz
- CE Class: IIa (directive 93/42/EEC)

### ADDITIONAL INFORMATION

Animal	Mixed
SKU	129-600-accessories
Overall Product Dimensions (mm)	370H x 210W x 350L

**Generic Size**

Medium

**Colour/Style**

White

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights