

NON VACUUM LONG CHAMBER AUTOCLAVE

2560-V S3 A



PRODUCT DESCRIPTION

The Vetclave 2560 Long Chamber Model has been designed specifically to meet the needs of your practice.

You have the ability to increase the sterilisation and drying times to suit the load. A very important feature for ensuring all bags and more complex instruments are thoroughly sterilised.

The large chamber size of 30-litre capacity enables larger instruments to be sterilised, this is ideal for equine and large animal instruments. The Instaclave has space for 6 trays. 3 trays are included as well as your hose, drip tray and handle for ease of use. The door gasket can be replaced and the overheat cut out can be reset easily by the user, reducing your maintenance costs and downtime.

Why choose a non-vacuum Autoclave?

A non-vacuum autoclave is designed specifically for sterilisation of non-wrapped and solid instruments. It comes with the added benefit of being able to increase the sterilising and drying times to suit the load.

Integrated Safety Features

You will always feel comfortable whilst using the Instaclave due to its abundance of safety features. The manual and mechanical door locking mechanism ensures that it cannot be opened whilst sterilising due to the pressure forcing the second lock to engage. Overheating protection and a pressure safety valve are two other features that will guarantee you an easy and safe experience.

Key Features:

- Ability to increase sterilising and drying times

- 2 Cycle types
- Larger chamber
- Safety features including manual and mechanical door locks and overheating protection
- Easy to use, clear buttons to switch between cycles, with clear lights to show your selection
- Non recycled water system

Technical Specifications	
Cycle Types	134°C Non Vacuum Cycle (drying cycle available) 121°C Non Vacuum Cycle (drying cycle available)
Chamber Size	30 litres
Capacity	6 Trays (Supplied With 3)
Dimensions	2560-V 72cm x 54cm x 41cm DWH
Weight	53kg
Cycle Times	134°C Non-Vacuum Cycle (30 or 40mins depending on drying cycle)* 121°C Non-Vacuum Cycle (30 or 40mins depending on drying cycle)*
Power Requirements	220/240v, Current: 9 amps

ADDITIONAL INFORMATION

SKU	2560-V S3 A
Overall Product Dimensions (mm)	720 x 540 x 410 DWH
Colour/Style	Stainless Steel
Brand	Burtons
Number of Cycles	2
Vacuum	No
Chamber Size	30 Litres
Autoclave Accessory/Spares	Yes

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights

