

## U500 - HEATED ULTRAWAVE CLEANING BATH - 4.45L

116-606



### PRODUCT DESCRIPTION

## The U500 Heated ultrasonic cleaning bath is ideal for cleaning small items in

#### Includes a FREE Pack of 6 Ultraclean M2 Pack Detergent

A heated ultrasonic bath with an analog timer and temperature control, the U500H has a 4.5 capacity tank. The heated version will speed up the ultrasonic cleaning process. Simply fill the tank with water, turn the dial to the required minutes, temperature and hit the Sonics button.

The U500H is supplied complete with a stainless steel basket and an ABS plastic lid. The lid has been designed to be used as a drip tray to collect excess water when the basket is removed

The U500H provides the cost benefits of excellent entry level ultrasonics with Ultrawave's long renowned quality and reliability.

Please note: 3 year warranty included when you order with us\*

*\*1 year as standard, 2 years extra warranty applicable if registered online with Ultrawave within 30 days of placing your order*

|                          |                             |
|--------------------------|-----------------------------|
| Overall Dimensions       | 330 x 180 x 275 (W x D x H) |
| Internal Tank Dimensions | 300 x 150 x 150mm           |
| Basket Dimensions        | 270 x 125 x 90mm            |

|                     |  |
|---------------------|--|
| Capacity            | Maximum: 5 Litres   Working: 4.5 Litres                                  |
| Construction        | 304 Stainless Steel Tank   |
| Control             | Digital Control<br>Time (0 - 15minutes)<br>Temperature (Ambient to 70°C) |
| Ultrasonic Power    | 100W   |
| Operating Frequency | 32.5-37.5KHz   |
| Heating Power       | 150W   |
| Required Services   | 13amp 3 Pin Socket   |

## ADDITIONAL INFORMATION

|                                 |           |
|---------------------------------|-----------|
| <b>Animal</b>                   | N/A       |
| <b>SKU</b>                      | 116-606   |
| <b>Brand</b>                    | Ultrawave |
| <b>Ultrasonic Bath Controls</b> | Analogue  |
| <b>Ultrasonic Bath Litres</b>   | 4.45      |

### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights