

**ECONOMY ULTRASONIC BATH 4.5L**

116-650

**PRODUCT DESCRIPTION**

The economy Ultrasonic Bath has a 5ltr capacity tank and is ideal for cleaning small items

Overall Dimensions	550 x 260 x 340mm (W x D x H)
Internal Tank Dimensions	240 x 142 x 149mm (W x D x H)
Basket Dimensions	210 x 120 x 110mm (W x D x H)
Capacity	5 Litres
Control	Digital Control
Ultrasonic Power	120W
Operating Frequency	40kHz
Required Services	13amp 3 Pin Socket

**ADDITIONAL INFORMATION**

<b>Animal</b>	Mixed
<b>SKU</b>	116-650
<b>Overall Product Dimensions (mm)</b>	550 x 260 x 340mm (W x D x H)

## Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

**Price Matching:** We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

**Backorders:** All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

**Returns Policy:** You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights