

WATER DISTILLER 4L WHITE STAINLESS STEEL WITH POLY JUG

139-200

**PRODUCT DESCRIPTION**

This automatic water distiller is very easy to use, produces over 4 litres of distilled water in only a few hours, and can be left running overnight, ensuring a full supply of distilled water for use with your autoclave each day

This unit can produce over 16 litres of distilled water every 24 hours at an average cost of 8p per litre.

No plumbing is required

Just fill with tap water and turn on the power. The water produced is ideal for use with your autoclave, dental system or other equipment requiring pure, distilled water. Why store large quantities of water in bulky, heavy containers when a distiller will provide you with consistently pure, clean water every day, taking up a fraction of your valuable space.

Water Distiller Product Specification

Water Distiller Output:	16 Litres per day
Size:	20cm wide x 36cm tall
Power source:	240V
Power consumption	580W (total)

Heater:	565W/Cooling Fan 15W
Weight:	3.5Kg (7.7lbs)
Capacity:	4lt per 3-4 hours

ADDITIONAL INFORMATION

Animal	Mixed
SKU	139-200

Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights