

LABCOLD TEMPERATURE LOGGER

125-046

**PRODUCT DESCRIPTION****Convenient, user-friendly USB data logger with wide application that works just like a removable memory drive**

Small and light, the logger plugs directly into a USB port of any suitable computer for configuration, data download and analysis. Up to 32,000 data points can be recorded before download is necessary.

ADDITIONAL INFORMATION

| | |
|---------------------|---------|
| Animal | N/A |
| SKU | 125-046 |
| Colour/Style | Orange |
| Brand | Labcold |

Good to know...

At Burton's we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs.

Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it!

Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates.

Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights