

## LABCOLD TEMPERATURE LOGGER

125-046



### PRODUCT DESCRIPTION

---

## Convenient, alarmed, user-friendly USB 2.0 data logger with wide application that works just like a removable memory drive

These loggers have been specially designed for the monitoring of the pharmaceutical cold chain and are ideal for using when vaccines or in transit or as an independent thermometer for your vaccine refrigerator. Labcold data loggers are specially designed to help with cold change compliance.

Small, lightweight, and portable, these waterproof loggers plug directly into a USB port of any suitable computer for configuration, data download and analysis using free software supplied with the logger. Up to 32,000 data points can be recorded before download is necessary. The report generated by the software also includes a data curve, a data report and every single measurement which can be stored to your PC or printed as a hard copy.

Easy to use, you can simply start the logger using the preprogrammed parameters described on the reverse of the logger or fully configure, including applying password protection by using the free downloadable software.

The software allows you to configure up to 5 separate alarm parameters and the PDF report generated by the loggers will colour code every time alarm limits are exceeded. If temperature deviation happens, it alarms you, with an LED screen which can display alarm alerts in real time which minimises the risks during the monitoring process.

Simply put the temperature loggers in with the product and once initiated they will record the temperature automatically with no more intervention. Our transport loggers also have a handy light that will flash green if the product stays within range, but flash red should there be a problem. These loggers can be used time and time again as they have an easy change battery.

The report generated by the software also includes a data curve, a data report and every single measurement which can be stored to your PC or printed as a hard copy.

The logger comes with a 6-point manufacturers' validation certificate but for extra peace of mind you can have it UKAS traceable calibrated to 3 points.

### Specification

- 2.0 A-Type USB connection
- Accuracy: +/-0.5 at -20°C to +40°C
- Up to 5 alarm settings programmable
- Adjustable sample rate
- Replaceable battery
- Visual alarm (if configured)
- Creates a PDF report which includes a data curve, data report and all data points
- FREE software - Windows compatible only
- Capacity: 32,000 data points
- 80mm (H) x 34mm (W) x 14mm (D)
- Weight: 25g
- Temperature Range: -30°C to 70°C

### ADDITIONAL INFORMATION

<b>SKU</b>	125-046
<b>Colour/Style</b>	Orange
<b>Brand</b>	Labcold

#### Good to know...

At Burtons we know you have a choice to shop with us, that is why we always aim to delight our customers with a fast, easy, shopping experience. Making and sourcing only the best products we can find to fit your needs. Price Matching: We regularly check all our prices against competitors so you don't have to. If you do find a like for like product cheaper we will aim to not just match it, but beat it! Backorders: All back ordered items will be ordered and shipped as quickly as we can to you. Special order items may take up to 4 weeks. Once ordered a member of our team will be in contact to keep you informed of estimated delivery dates. Returns Policy: You may return new, unworn or unused products within fourteen (14) days of delivery for a full refund of the cost of the goods, or an exchange if preferred. This is on top of your statutory rights